

Resource Summary Report

Generated by [RRID](#) on Sep 8, 2026

ZIC2 antibody

RRID:AB_371699

Type: Antibody

Proper Citation

GeneTex Cat# GTX12072, RRID:AB_371699

Antibody Information

URL: http://antibodyregistry.org/AB_371699

Proper Citation: GeneTex Cat# GTX12072, RRID:AB_371699

Target Antigen: ZIC2 antibody

Host Organism: rabbit

Clonality: polyclonal

Comments: Discontinued; manufacturer recommendations: IgG; IgG IHC: Use at a dilution of 1/200 - 1/1000. WB: Use at a dilution of 1/200 - 1/1000. Detects a band of approximately 55 kDa. Preparation Instructions Method 1 for immunostaining and immunoblotting (Western blot)1. Carefully resuspend antibody pellet to uniformity.2. Remove a fixed amount of suspension and dissolve 1/5 in PBS or TBS to yield a 100µg/ml solution. Method 2 for immunoprecipitation, supershift, immunostaining and immunoblotting (Western blot)1 Pellet antibodies at 10,000 to 15,000 x g for 10 minutes at 2 to 8 °C using a microcentrifuge.2 Carefully remove as much supernatant as possible. It is not necessary to remove all the ammonium sulfate solution; a small residual amount will not effect the antibody preparation. Dissolve the pellet (antibody) in small volume (50µ

Antibody Name: ZIC2 antibody

Description: This polyclonal targets ZIC2 antibody

Target Organism: rat, mouse, human

Antibody ID: AB_371699

Vendor: GeneTex

Catalog Number: GTX12072

Record Creation Time: 20231110T081209+0000

Record Last Update: 20260830T032200+0000

Ratings and Alerts

No rating or validation information has been found for ZIC2 antibody.

Warning: Discontinued at GeneTex

Discontinued; manufacturer recommendations: IgG; IgG IHC: Use at a dilution of 1/200 - 1/1000. WB: Use at a dilution of 1/200 - 1/1000. Detects a band of approximately 55 kDa. Preparation Instructions Method 1 for immunostaining and immunoblotting (Western blot) 1. Carefully resuspend antibody pellet to uniformity. 2. Remove a fixed amount of suspension and dissolve 1/5 in PBS or TBS to yield a 100µg/ml solution. Method 2 for immunoprecipitation, supershift, immunostaining and immunoblotting (Western blot) 1 Pellet antibodies at 10,000 to 15,000 x g for 10 minutes at 2 to 8 °C using a microcentrifuge. 2 Carefully remove as much supernatant as possible. It is not necessary to remove all the ammonium sulfate solution; a small residual amount will not effect the antibody preparation. Dissolve the pellet (antibody) in small volume (50µ

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.