

Resource Summary Report

Generated by [RRID](#) on Aug 11, 2026

Filamin [FLMN01 (PM6/317)] antibody

RRID:AB_384986

Type: Antibody

Proper Citation

GeneTex Cat# GTX23261, RRID:AB_384986

Antibody Information

URL: http://antibodyregistry.org/AB_384986

Proper Citation: GeneTex Cat# GTX23261, RRID:AB_384986

Target Antigen: Filamin [FLMN01 (PM6/317)] antibody

Host Organism: mouse

Clonality: monoclonal

Comments: Discontinued; manufacturer recommendations: Immunohistochemistry - fixed; Western Blot; Immunohistochemistry; IHC-P, WB, IHC (Formalin-fixed paraffin-embedded sections), Western blot. The usefulness of this product in other applications has not been determined., Western Blotting: use at a concentration of 1 - 2µg / ml for 2hrs at RT. Immunohistochemistry (formalin fixed paraffin embedded): use at a concentration of 1 - 2µg / ml for 30 min at RT. Not tested in other applications. Optimal dilutions / concentrations should be determined by the end user. By Western blot, this antibody detects a band of 250 kDa, which corresponds to the predicted molecular weight of Filamin. Immunohistochemical staining of formalin / paraffin tissues requires boiling tissue sections in 10mM citrate buffer, pH 6.0, for 10 - 20 min followed by

Antibody Name: Filamin [FLMN01 (PM6/317)] antibody

Description: This monoclonal targets Filamin [FLMN01 (PM6/317)] antibody

Target Organism: rat, mouse, chickenbird, rabbit, human

Antibody ID: AB_384986

Vendor: GeneTex

Catalog Number: GTX23261

Record Creation Time: 20250819T205728+0000

Ratings and Alerts

No rating or validation information has been found for Filamin [FLMN01 (PM6/317)] antibody.

Warning: Discontinued at GeneTex

Discontinued; manufacturer recommendations: Immunohistochemistry - fixed; Western Blot; Immunohistochemistry; IHC-P, WB, IHC (Formalin-fixed paraffin-embedded sections), Western blot. The usefulness of this product in other applications has not been determined., Western Blotting: use at a concentration of 1 - 2µg / ml for 2hrs at RT. Immunohistochemistry (formalin fixed paraffin embedded): use at a concentration of 1 - 2µg / ml for 30 min at RT. Not tested in other applications. Optimal dilutions / concentrations should be determined by the end user. By Western blot, this antibody detects a band of 250 kDa, which corresponds to the predicted molecular weight of Filamin. Immunohistochemical staining of formalin / paraffin tissues requires boiling tissue sections in 10mM citrate buffer, pH 6.0, for 10 - 20 min followed by

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.