

Resource Summary Report

Generated by [RRID](#) on Sep 10, 2026

CD15 [MMA] antibody (ready-to-use)

RRID:AB_367200

Type: Antibody

Proper Citation

GeneTex Cat# GTX27495, RRID:AB_367200

Antibody Information

URL: http://antibodyregistry.org/AB_367200

Proper Citation: GeneTex Cat# GTX27495, RRID:AB_367200

Target Antigen: CD15 [MMA] antibody (ready-to-use)

Host Organism: mouse

Clonality: monoclonal

Comments: Discontinued; manufacturer recommendations: Immunohistochemistry - frozen; Immunohistochemistry; Immunohistochemistry - fixed; IHC-P: Use at an assay dependent dilution. IHC-Fr: Use at an assay dependent dilution. We suggest an incubation period of 30 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the end user. Formalin fixed paraffin embedded tissue sections require high temperature antigen unmasking with 10 mM citrate buffer, pH 6.0 prior to immunostaining. Not tested in other applications. Optimal dilutions/concentrations should be determined by the end user., IHC (Frozen sections), IHC (Formalin-fixed paraffin-embedded sections). The usefulness of this product in other applications has not been determined., IHC-Fr, IHC-P

Antibody Name: CD15 [MMA] antibody (ready-to-use)

Description: This monoclonal targets CD15 [MMA] antibody (ready-to-use)

Target Organism: mouse, human - not yet tested in other species, human

Antibody ID: AB_367200

Vendor: GeneTex

Catalog Number: GTX27495

Record Creation Time: 20231110T081250+0000

Ratings and Alerts

No rating or validation information has been found for CD15 [MMA] antibody (ready-to-use).

Warning: Discontinued at GeneTex

Discontinued; manufacturer recommendations: Immunohistochemistry - frozen; Immunohistochemistry; Immunohistochemistry - fixed; IHC-P: Use at an assay dependent dilution. IHC-Fr: Use at an assay dependent dilution. We suggest an incubation period of 30 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the end user. Formalin fixed paraffin embedded tissue sections require high temperature antigen unmasking with 10 mM citrate buffer, pH 6.0 prior to immunostaining. Not tested in other applications. Optimal dilutions/concentrations should be determined by the end user., IHC (Frozen sections), IHC (Formalin-fixed paraffin-embedded sections). The usefulness of this product in other applications has not been determined., IHC-Fr, IHC-P

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.