

# Resource Summary Report

Generated by [RRID](#) on Aug 5, 2026

## GRB2 [2GB04] antibody

RRID:AB\_368127

Type: Antibody

### Proper Citation

GeneTex Cat# GTX23655, RRID:AB\_368127

### Antibody Information

**URL:** [http://antibodyregistry.org/AB\\_368127](http://antibodyregistry.org/AB_368127)

**Proper Citation:** GeneTex Cat# GTX23655, RRID:AB\_368127

**Target Antigen:** GRB2 [2GB04] antibody

**Host Organism:** mouse

**Clonality:** monoclonal

**Comments:** Discontinued; manufacturer recommendations: Immunohistochemistry - fixed; Immunohistochemistry - frozen; Immunohistochemistry; Immunohistochemistry (formalin fixed paraffin embedded / frozen sections): use at a dilution of 1/50 - 1/75 in an ABC method. Not tested in other applications. Optimal dilutions / concentrations should be determined by the end user. Formalin fixed paraffin embedded tissue sections require high temperature antigen unmasking with 10 mM citrate buffer, pH 6.0 prior to immunostaining., IHC (Frozen sections), IHC (Formalin-fixed paraffin-embedded sections). The usefulness of this product in other applications has not been determined., IHC-Fr, IHC-P

**Antibody Name:** GRB2 [2GB04] antibody

**Description:** This monoclonal targets GRB2 [2GB04] antibody

**Target Organism:** rat, xenopusamphibian, mouse, bovine, human

**Antibody ID:** AB\_368127

**Vendor:** GeneTex

**Catalog Number:** GTX23655

**Record Creation Time:** 20250820T033000+0000

**Record Last Update:** 20250820T162449+0000

---

## Ratings and Alerts

No rating or validation information has been found for GRB2 [2GB04] antibody.

**Warning:** Discontinued at GeneTex

Discontinued; manufacturer recommendations: Immunohistochemistry - fixed; Immunohistochemistry - frozen; Immunohistochemistry; Immunohistochemistry (formalin fixed paraffin embedded / frozen sections): use at a dliution of 1/50 - 1/75 in an ABC method. Not tested in other applications. Optimal dilutions / concentrations should be determined by the end user. Formalin fixed paraffin embedded tissue sections require high temperature antigen unmasking with 10 mM citrate buffer, pH 6.0 prior to immunostaining., IHC (Frozen sections), IHC (Formalin-fixed paraffin-embedded sections). The usefulness of this product in other applications has not been determined., IHC-Fr, IHC-P

---

## Data and Source Information

**Source:** [Antibody Registry](#)

---

## Usage and Citation Metrics

We have not found any literature mentions for this resource.