

Resource Summary Report

Generated by [RRID](#) on Sep 11, 2026

Melanoma [HMB45] antibody (ready-to-use)

RRID:AB_369646

Type: Antibody

Proper Citation

GeneTex Cat# GTX20906, RRID:AB_369646

Antibody Information

URL: http://antibodyregistry.org/AB_369646

Proper Citation: GeneTex Cat# GTX20906, RRID:AB_369646

Target Antigen: Melanoma [HMB45] antibody (ready-to-use)

Host Organism: mouse

Clonality: monoclonal

Comments: Discontinued; manufacturer recommendations: Immunohistochemistry - frozen; Immunohistochemistry - fixed; Immunohistochemistry; This antibody has been pretitered and quality controlled to work allon formalin-fixed paraffin-embedded as well as acetone fixed cryostat tissue sections. No further titration is required. We suggest an incubation period of 30-60 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the user. Enzymatic predigestion with proteolytic enzymes may increase staining intensity., IHC (Frozen sections), IHC (Formalin-fixed paraffin-embedded sections). The usefulness of this product in other applications has not been determined., IHC-Fr, IHC-P

Antibody Name: Melanoma [HMB45] antibody (ready-to-use)

Description: This monoclonal targets Melanoma [HMB45] antibody (ready-to-use)

Antibody ID: AB_369646

Vendor: GeneTex

Catalog Number: GTX20906

Record Creation Time: 20231110T081254+0000

Record Last Update: 20260829T002348+0000

Ratings and Alerts

No rating or validation information has been found for Melanoma [HMB45] antibody (ready-to-use).

Warning: Discontinued at GeneTex

Discontinued; manufacturer recommendations: Immunohistochemistry - frozen; Immunohistochemistry - fixed; Immunohistochemistry; This antibody has been pretitered and quality controlled to work allon formalin-fixed paraffin-embedded as well as acetone fixed cryostat tissue sections. No further titration is required. We suggest an incubation period of 30-60 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the user. Enzymatic predigestion with proteolytic enzymes may increase staining intensity., IHC (Frozen sections), IHC (Formalin-fixed paraffin-embedded sections). The usefulness of this product in other applications has not been determined., IHC-Fr, IHC-P

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.