

Resource Summary Report

Generated by [RRID](#) on Aug 10, 2026

Proteosome Activator, 200 kDa antibody

RRID:AB_372005

Type: Antibody

Proper Citation

GeneTex Cat# GTX25620, RRID:AB_372005

Antibody Information

URL: http://antibodyregistry.org/AB_372005

Proper Citation: GeneTex Cat# GTX25620, RRID:AB_372005

Target Antigen: Proteosome Activator 200 kDa antibody

Host Organism: rabbit

Clonality: polyclonal

Comments: Discontinued; manufacturer recommendations: IgG; IgG Immunocytochemistry, Immunofluorescence, Immunoprecipitation, Western Blot. IF: Use at a dilution of 1/600. For ICC/IF staining of PA200 in HeLa cells results in diffuse nuclear staining prior to ionizing radiation treatment. After ionizing radiation treatment it results in focal aggregate staining within the nucleus. For IP: Use at an assay dependent dilution. WB: Use at a dilution of 1/1000. Detects a band of approximately 200 kDa representing PA200 from bovine testis tissue extract. In mouse tissues, this antibody detects both a 200 kDa and a 160 kDa protein. In mouse tissues, the 200 kDa species is most abundant in testis whereas the 160 kDa species is enriched in the spleen. Not tested in other applications. Optimal dilutions/concentrations should be determined by the end user.; Im

Antibody Name: Proteosome Activator, 200 kDa antibody

Description: This polyclonal targets Proteosome Activator 200 kDa antibody

Target Organism: rat, mouse, bovine, human

Antibody ID: AB_372005

Vendor: GeneTex

Catalog Number: GTX25620

Record Creation Time: 20250820T011217+0000

Ratings and Alerts

No rating or validation information has been found for Proteasome Activator, 200 kDa antibody.

Warning: Discontinued at GeneTex

Discontinued; manufacturer recommendations: IgG; IgG Immunocytochemistry, Immunofluorescence, Immunoprecipitation, Western Blot. IF: Use at a dilution of 1/600. For ICC/IF staining of PA200 in HeLa cells results in diffuse nuclear staining prior to ionizing radiation treatment. After ionizing radiation treatment it results in focal aggregate staining within the nucleus. For IP: Use at an assay dependent dilution. WB: Use at a dilution of 1/1000. Detects a band of approximately 200 kDa representing PA200 from bovine testis tissue extract. In mouse tissues, this antibody detects both a 200 kDa and a 160 kDa protein. In mouse tissues, the 200 kDa species is most abundant in testis whereas the 160 kDa species is enriched in the spleen. Not tested in other applications. Optimal dilutions/concentrations should be determined by the end user.; Im

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.