

Resource Summary Report

Generated by [RRID](#) on Jul 29, 2026

VEGF antibody

RRID:AB_373277

Type: Antibody

Proper Citation

GeneTex Cat# GTX22992, RRID:AB_373277

Antibody Information

URL: http://antibodyregistry.org/AB_373277

Proper Citation: GeneTex Cat# GTX22992, RRID:AB_373277

Target Antigen: VEGF antibody

Host Organism: rabbit

Clonality: polyclonal

Comments: Discontinued; manufacturer recommendations: IgG; IgG Immunohistochemistry, IHC (Formalin-fixed paraffin-embedded sections). Immunohistochemistry (formalin fixed paraffin embedded, rodent tissues only): use at a concentration of 5 - 10µg / ml for 30 min at RT. Staining of formalin fixed tissues requires boiling tissue sections in 10mM Tris-HCl, pH 10.0, for 10 - 20 min followed by cooling at RT for 20 min. Neutralising (Int J Oncology, 8:505-511): assay dependent. Not tested in other applications. Optimal dilutions / concentrations should be determined by the end user. Staining of formalin fixed tissues requires boiling tissue sections in 10mM Tris-HCl, pH 10.0, for 10 - 20 min followed by cooling at RT for 20 min.; Block/Neutralize/Inhibit; Immunohistochemistry; Immunohistochemistry - fixed

Antibody Name: VEGF antibody

Description: This polyclonal targets VEGF antibody

Target Organism: rat, mouse, human

Antibody ID: AB_373277

Vendor: GeneTex

Catalog Number: GTX22992

Record Creation Time: 20250819T205844+0000

Ratings and Alerts

No rating or validation information has been found for VEGF antibody.

Warning: Discontinued at GeneTex

Discontinued; manufacturer recommendations: IgG; IgG Immunohistochemistry, IHC (Formalin-fixed paraffin-embedded sections). Immunohistochemistry (formalin fixed paraffin embedded, rodent tissues only): use at a concentration of 5 - 10µg / ml for 30 min at RT. Staining of formalin fixed tissues requires boiling tissue sections in 10mM Tris-HCl, pH 10.0, for 10 - 20 min followed by cooling at RT for 20 min. Neutralising (Int J Oncology, 8:505-511): assay dependent. Not tested in other applications. Optimal dilutions / concentrations should be determined by the end user. Staining of formalin fixed tissues requires boiling tissue sections in 10mM Tris-HCl, pH 10.0, for 10 - 20 min followed by cooling at RT for 20 min.; Block/Neutralize/Inhibit; Immunohistochemistry; Immunohistochemistry - fixed

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.