

Resource Summary Report

Generated by [RRID](#) on Sep 10, 2026

Cardiac Troponin I [2D5] antibody

RRID:AB_374201

Type: Antibody

Proper Citation

GeneTex Cat# GTX72420, RRID:AB_374201

Antibody Information

URL: http://antibodyregistry.org/AB_374201

Proper Citation: GeneTex Cat# GTX72420, RRID:AB_374201

Target Antigen: Cardiac Troponin I [2D5] antibody

Host Organism: mouse

Clonality: monoclonal

Comments: Discontinued; manufacturer recommendations: Immunohistochemistry; Immunohistochemistry - frozen; Other; Immunohistochemistry - fixed; Immunohistochemistry 1:25-1:50.(For formalin-fixed paraffin-embedded (FFPE) tissue sections antigen demasker (10 mM citrate buffer pH 6.0) is required. Boil the tissue 10-20 mins in this buffer, followed by cooling to RT. This antibody can be used on formalin-fixed, paraffin-embedded tissue sections. Prolonged fixation in buffered formalin can destroy the epitope. The antibody may be used at a dilution of 1:25-1:50 with LAB/Probe. When using paraffin-embedded and formalin-fixed tissues, high temperature antigenic unmasking technique (see instructions below) is strongly recommended for consistent and reproducible results. This product can also be used on frozen tissue sections. The optimal conditions

Antibody Name: Cardiac Troponin I [2D5] antibody

Description: This monoclonal targets Cardiac Troponin I [2D5] antibody

Antibody ID: AB_374201

Vendor: GeneTex

Catalog Number: GTX72420

Record Creation Time: 20231110T081215+0000

Record Last Update: 20260831T213452+0000

Ratings and Alerts

No rating or validation information has been found for Cardiac Troponin I [2D5] antibody.

Warning: Discontinued at GeneTex

Discontinued; manufacturer recommendations: Immunohistochemistry;

Immunohistochemistry - frozen; Other; Immunohistochemistry - fixed;

Immunohistochemistry 1:25-1:50.(For formalin-fixed paraffin-embedded (FFPE) tissue sections antigen demasker (10 mM citrate buffer pH 6.0) is required. Boil the tissue 10-20 mins in this buffer, followed by cooling to RT. This antibody can be used on formalin-fixed, paraffin-embedded tissue sections. Prolonged fixation in buffered formalin can destroy the epitope. The antibody may be used at a dilution of 1:25-1:50 with LAB/Probe. When using paraffin-embedded and formalin-fixed tissues, high temperature antigenic unmasking technique (see instructions below) is strongly recommended for consistent and reproducible results. This product can also be used on frozen tissue sections. The optimal conditions

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.