

Resource Summary Report

Generated by [RRID](#) on Jul 29, 2026

CD68 [KP1] antibody (ready-to-use)

RRID:AB_374277

Type: Antibody

Proper Citation

GeneTex Cat# GTX20845, RRID:AB_374277

Antibody Information

URL: http://antibodyregistry.org/AB_374277

Proper Citation: GeneTex Cat# GTX20845, RRID:AB_374277

Target Antigen: CD68 [KP1] antibody (ready-to-use)

Host Organism: mouse

Clonality: monoclonal

Comments: Discontinued; manufacturer recommendations: This antibody has been pretitered and quality controlled to work allon formalin-fixed paraffin-embedded as well as acetone fixed cryostat tissue sections. No further titration is required. We suggest an incubation period of 60 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the user. Proteolytic treatment is required prior to immunostaining of formalin fixed paraffin embedded tissue sections. High temperature antigen unmasking with 10mM citrate buffer, pH6.0 prior to staining., IHC (Frozen sections), IHC (Formalin-fixed paraffin-embedded sections). The usefulness of this product in other applications has not been determined., IHC-Fr, IHC-P; Immunohistochemistry - fixed; Immunohistochem

Antibody Name: CD68 [KP1] antibody (ready-to-use)

Description: This monoclonal targets CD68 [KP1] antibody (ready-to-use)

Target Organism: human

Antibody ID: AB_374277

Vendor: GeneTex

Catalog Number: GTX20845

Record Creation Time: 20250820T022253+0000

Ratings and Alerts

No rating or validation information has been found for CD68 [KP1] antibody (ready-to-use).

Warning: Discontinued at GeneTex

Discontinued; manufacturer recommendations: This antibody has been pretitered and quality controlled to work allon formalin-fixed paraffin-embedded as well as acetone fixed cryostat tissue sections. No further titration is required. We suggest an incubation period of 60 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the user. Proteolytic treatment is required prior to immunostaining of formalin fixed paraffin embedded tissue sections. High temperature antigen unmasking with 10mM citrate buffer, pH6.0 prior to staining., IHC (Frozen sections), IHC (Formalin-fixed paraffin-embedded sections). The usefulness of this product in other applications has not been determined., IHC-Fr, IHC-P; Immunohistochemistry - fixed; Immunohistochem

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.