

Resource Summary Report

Generated by [RRID](#) on Aug 10, 2026

Proinsulin [1D3] antibody

RRID:AB_374546

Type: Antibody

Proper Citation

GeneTex Cat# GTX28300, RRID:AB_374546

Antibody Information

URL: http://antibodyregistry.org/AB_374546

Proper Citation: GeneTex Cat# GTX28300, RRID:AB_374546

Target Antigen: Proinsulin [1D3] antibody

Host Organism: mouse

Clonality: monoclonal

Comments: Discontinued; manufacturer recommendations: Radioimmunoassay; ELISA; ELISA, RIA, ELISA, Radioimmunoassay. The usefulness of this product in other applications has not been determined., This antibody is used for direct capture of prohormone from serum sample with insulin-specific GTX28302 (7F8) as a conjugate. This antibody and GTX28299 (1D6) can be also used for construction of a sensitive two-site sandwich enzyme-linked/RIA immunosorbent assays for the measurement of intact human proinsulin in serum or plasma with the detection limit in buffer 0.4 - 0.5 pmol/L. Specific immunosorbent (immobilized 1D3) removal of proinsulin molecules from serum samples before the measurement of insulin C-peptide.

Antibody Name: Proinsulin [1D3] antibody

Description: This monoclonal targets Proinsulin [1D3] antibody

Target Organism: human

Antibody ID: AB_374546

Vendor: GeneTex

Catalog Number: GTX28300

Record Creation Time: 20250819T215952+0000

Record Last Update: 20250820T004622+0000

Ratings and Alerts

No rating or validation information has been found for Proinsulin [1D3] antibody.

Warning: Discontinued at GeneTex

Discontinued; manufacturer recommendations: Radioimmunoassay; ELISA; ELISA, RIA, ELISA, Radioimmunoassay. The usefulness of this product in other applications has not been determined., This antibody is used for direct capture of prohormone from serum sample with insulin-specific GTX28302 (7F8) as a conjugate. This antibody and GTX28299 (1D6) can be also used for construction of a sensitive two-site sandwich enzyme-linked/RIA immunosorbent assays for the measurement of intact human proinsulin in serum or plasma with the detection limit in buffer 0.4 - 0.5 pmol/L. Specific immunosorbent (immobilized 1D3) removal of proinsulin molecules from serum samples before the measurement of insulinic C-peptide.

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.