

Resource Summary Report

Generated by [RRID](#) on Aug 11, 2026

Ki67 antibody

RRID:AB_375075

Type: Antibody

Proper Citation

GeneTex Cat# GTX73691, RRID:AB_375075

Antibody Information

URL: http://antibodyregistry.org/AB_375075

Proper Citation: GeneTex Cat# GTX73691, RRID:AB_375075

Target Antigen: Ki67 antibody

Host Organism: rabbit

Clonality: polyclonal

Comments: Discontinued; manufacturer recommendations: IgG; IgG Immunohistochemistry, IHC (Formalin-fixed paraffin-embedded sections). This antibody may be diluted to a titer of 1:25-1:50 in an ABC method. However, depending on the staining conditions employed, we suggest that final dilution should be determined by the user. We suggest an incubation period of 30 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation conditions should be determined by the user. Formalin fixed paraffin embedded tissue sections require high temperature antigen unmasking with 10 mM citrate buffer, pH 6.0 prior to immunostaining. Not tested; Immunohistochemistry; Immunohistochemistry - fixed

Antibody Name: Ki67 antibody

Description: This polyclonal targets Ki67 antibody

Target Organism: feline, rat, hamster, porcine, donkey, canine, goat, horse, mouse, rabbit, bovine, human, sheep

Antibody ID: AB_375075

Vendor: GeneTex

Catalog Number: GTX73691

Record Creation Time: 20250820T052716+0000

Record Last Update: 20250820T213640+0000

Ratings and Alerts

No rating or validation information has been found for Ki67 antibody.

Warning: Discontinued at GeneTex

Discontinued; manufacturer recommendations: IgG; IgG Immunohistochemistry, IHC (Formalin-fixed paraffin-embedded sections). This antibody may be diluted to a titer of 1:25-1:50 in an ABC method. However, depending on the staining conditions employed, we suggest that final dilution should be determined by the user. We suggest an incubation period of 30 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation conditions should be determined by the user. Formalin fixed paraffin embedded tissue sections require high temperature antigen unmasking with 10 mM citrate buffer, pH 6.0 prior to immunostaining. Not tested; Immunohistochemistry; Immunohistochemistry - fixed

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.