

Resource Summary Report

Generated by RRID on Aug 6, 2026

Horiba: Camsizer X2 Dynamic Image Analysis System Particle Size

RRID:SCR_020040

Type: Tool

Proper Citation

Horiba: Camsizer X2 Dynamic Image Analysis System Particle Size (RRID:SCR_020040)

Resource Information

URL: https://www.horiba.com/en_en/products/detail/action/show/Product/camsizer-x2-1133/

Proper Citation: Horiba: Camsizer X2 Dynamic Image Analysis System Particle Size (RRID:SCR_020040)

Description: Dynamic image analysis system that provides particle size and shape distribution analysis needed for quality control of fine powders and suspensions, measuring 0.8 μm to 8 mm.

Resource Type: instrument resource

Keywords: Horiba, Dynamic Image Analysis System Particle Size, Instrument Equipment, USEDit,

Availability: Commercially available

Resource Name: Horiba: Camsizer X2 Dynamic Image Analysis System Particle Size

Resource ID: SCR_020040

Alternate IDs: Model_Number_Camsizer X2

Record Creation Time: 20220129T080348+0000

Record Last Update: 20260801T190650+0000

Ratings and Alerts

No rating or validation information has been found for Horiba: Camsizer X2 Dynamic Image Analysis System Particle Size.

No alerts have been found for Horiba: Camsizer X2 Dynamic Image Analysis System Particle Size.

Data and Source Information

Source: [SciCrunch Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.