

Resource Summary Report

Generated by [RRID](#) on Sep 9, 2026

MUNle007-A

RRID:CVCL_C860

Type: Cell Line

Proper Citation

RRID:CVCL_C860

Cell Line Information

URL: https://web.expasy.org/cellosaurus/CVCL_C860

Proper Citation: RRID:CVCL_C860

Description: Cell line MUNle007-A is a Embryonic stem cell with a species of origin Homo sapiens (Human)

Sex: Female

Defining Citation: [PMID:15955829](#), [PMID:17572666](#), [PMID:29248004](#), [PMID:30650618](#)

Comments: From: Masaryk University; Brno; Czech Republic.

Category: Embryonic stem cell

Organism: Homo sapiens (Human)

Name: MUNle007-A

Synonyms: CCTL 14, CCTL-14, CCTL14, Center for Cell Therapy Line 14

Cross References: hPSCreg:MUNle007-A, IGRhCellID:CCTL-14%20p38, IGRhCellID:CCTL-14%20p49, SKIP:SKIP001846, Wikidata:Q54907084

ID: CVCL_C860

Record Creation Time: 20220427T220827+0000

Record Last Update: 20260905T181222+0000

Ratings and Alerts

No rating or validation information has been found for MUNle007-A.

No alerts have been found for MUNle007-A.

Data and Source Information

Source: [CelloSaurus](#)

Usage and Citation Metrics

We found 6 mentions in open access literature.

Listed below are recent publications. The full list is available at [RRID](#) .

Portakal T, et al. (2025) Lipopolysaccharide induces retention of E-cadherin in the endoplasmic reticulum and promotes hybrid epithelial-to-mesenchymal transition of human embryonic stem cells-derived expandable lung epithelial cells. Inflammation research : official journal of the European Histamine Research Society ... [et al.], 74(1), 82.

Chochola V, et al. (2025) TYRAY-Functionalized Alginate Bioinks for 3D Bioprinting Support Stem Cell Culture and Endothelial Network Formation. ACS biomaterials science & engineering, 11(12), 7368.

Mandal G, et al. (2022) Matrix enrichment by black phosphorus improves ionization and reproducibility of mass spectrometry of intact cells, peptides, and amino acids. Scientific reports, 12(1), 1175.

Golunova A, et al. (2021) Direct and Indirect Biomimetic Peptide Modification of Alginate: Efficiency, Side Reactions, and Cell Response. International journal of molecular sciences, 22(11).

Va??hara P, et al. (2018) Intact Cell Mass Spectrometry as a Quality Control Tool for Revealing Minute Phenotypic Changes of Cultured Human Embryonic Stem Cells. Stem cells translational medicine, 7(1), 109.

Gudernova I, et al. (2017) One reporter for in-cell activity profiling of majority of protein kinase oncogenes. eLife, 6.