

Resource Summary Report

Generated by [RRID](#) on Sep 10, 2026

Streptococcus, group A carbohydrate antibody [0431] (FITC)

RRID:AB_11173417

Type: Antibody

Proper Citation

GeneTex Cat# GTX36638, RRID:AB_11173417

Antibody Information

URL: http://antibodyregistry.org/AB_11173417

Proper Citation: GeneTex Cat# GTX36638, RRID:AB_11173417

Target Antigen: Streptococcus group A carbohydrate antibody [0431] (FITC)

Host Organism: mouse

Clonality: monoclonal

Comments: Discontinued; manufacturer recommendations: Immunocytochemistry; Immunohistochemistry - fixed; Immunohistochemistry; Immunoprecipitation; Immunofluorescence; Recommended Starting Dilutions: For ICC/IF: Use at a starting dilution of 1:5 - 1:20. Acetone fixation is recommended.; Not yet tested in other applications. Optimal dilutions should be determined experimentally by the researcher.

Clone 0431 immunoprecipitates purified group A carbohydrate from both hot HCl & nitrous acid extracts. Works well HRP & gold conjugated. This monoclonal functions as a capture antibody on nitrocellulose but not on polystyrene solid phases. Clone 0721 immunoprecipitates purified group B carbohydrate from both autoclaved and enzyme prepared extracts but not HCl extracts., Immunoprecipitation. The usefulness of this product in other applications has

Antibody Name: Streptococcus, group A carbohydrate antibody [0431] (FITC)

Description: This monoclonal targets Streptococcus group A carbohydrate antibody [0431] (FITC)

Target Organism: bacteria

Antibody ID: AB_11173417

Vendor: GeneTex

Catalog Number: GTX36638

Record Creation Time: 20231110T060235+0000

Record Last Update: 20260829T045907+0000

Ratings and Alerts

No rating or validation information has been found for Streptococcus, group A carbohydrate antibody [0431] (FITC).

Warning: Discontinued at GeneTex

Discontinued; manufacturer recommendations: Immunocytochemistry; Immunohistochemistry - fixed; Immunohistochemistry; Immunoprecipitation; Immunofluorescence; Recommended Starting Dilutions: For ICC/IF: Use at a starting dilution of 1:5 - 1:20. Acetone fixation is recommended.; Not yet tested in other applications. Optimal dilutions should be determined experimentally by the researcher.

Clone 0431 immunoprecipitates purified group A carbohydrate from both hot HCl & nitrous acid extracts. Works well HRP & gold conjugated. This monoclonal functions as a capture antibody on nitrocellulose but not on polystyrene solid phases. Clone 0721 immunoprecipitates purified group B carbohydrate from both autoclaved and enzyme prepared extracts but not HCl extracts., Immunoprecipitation. The usefulness of this product in other applications has

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.