

Resource Summary Report

Generated by [RRID](#) on Jul 30, 2026

Donkey anti-Rat IgG (H+L) Cross-Adsorbed Secondary Antibody, DyLight™ 488

RRID:AB_2556606

Type: Antibody

Proper Citation

Thermo Fisher Scientific Cat# SA5-10026, RRID:AB_2556606

Antibody Information

URL: http://antibodyregistry.org/AB_2556606

Proper Citation: Thermo Fisher Scientific Cat# SA5-10026, RRID:AB_2556606

Target Antigen: Rat IgG (H+L)

Host Organism: donkey

Clonality: polyclonal secondary

Comments: Applications: ICC/IF (1:50-1:500), IHC (1:50-1:500), IP (Assay-dependent), WB (1:5,000-1:20,000), Flow (1:50-1:200)

Antibody Name: Donkey anti-Rat IgG (H+L) Cross-Adsorbed Secondary Antibody, DyLight™ 488

Description: This polyclonal secondary targets Rat IgG (H+L)

Target Organism: rat

Antibody ID: AB_2556606

Vendor: Thermo Fisher Scientific

Catalog Number: SA5-10026

Record Creation Time: 20250820T035327+0000

Record Last Update: 20250820T172723+0000

Ratings and Alerts

No rating or validation information has been found for Donkey anti-Rat IgG (H+L) Cross-Adsorbed Secondary Antibody, DyLight™ 488.

No alerts have been found for Donkey anti-Rat IgG (H+L) Cross-Adsorbed Secondary Antibody, DyLight™ 488.

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We found 7 mentions in open access literature.

Listed below are recent publications. The full list is available at [RRID](#) .

Puente-Ruiz SC, et al. (2025) B cell-derived nociceptin/orphanin FQ contributes to impaired glucose tolerance and insulin resistance in obesity. *iScience*, 28(7), 112819.

Nahiyera M, et al. (2025) Sca-1 expression depicts pro-inflammatory murine neutrophils under steady-state and pathological conditions. *Cell reports*, 44(7), 115995.

Zhang L, et al. (2024) Regulation of muscle hypertrophy through granulin: Relayed communication among mesenchymal progenitors, macrophages, and satellite cells. *Cell reports*, 43(4), 114052.

Sun F, et al. (2024) AdipoRon promotes amyloid-? clearance through enhancing autophagy via nuclear GAPDH-induced sirtuin 1 activation in Alzheimer's disease. *British journal of pharmacology*, 181(17), 3039.

Kaneshige A, et al. (2022) Relayed signaling between mesenchymal progenitors and muscle stem cells ensures adaptive stem cell response to increased mechanical load. *Cell stem cell*, 29(2), 265.

Winter C, et al. (2018) Chrono-pharmacological Targeting of the CCL2-CCR2 Axis Ameliorates Atherosclerosis. *Cell metabolism*, 28(1), 175.

Gabunia K, et al. (2016) IL-19 Halts Progression of Atherosclerotic Plaque, Polarizes, and Increases Cholesterol Uptake and Efflux in Macrophages. *The American journal of pathology*, 186(5), 1361.