

# Resource Summary Report

Generated by [RRID](#) on Aug 1, 2026

## Anti GAPDH in vitro Antibody

RRID:AB\_2858176

Type: Antibody

---

### Proper Citation

(Hytest Cat# 5G4cc-6C5cc, RRID:AB\_2858176)

---

### Antibody Information

**URL:** [http://antibodyregistry.org/AB\\_2858176](http://antibodyregistry.org/AB_2858176)

**Proper Citation:** (Hytest Cat# 5G4cc-6C5cc, RRID:AB\_2858176)

**Target Antigen:** GAPDH

**Host Organism:** mouse

**Clonality:** monoclonal

**Comments:** Applications: EIA, WB, IF, IHC, IP, porcine, canine, rabbit, cat, rat, mouse

**Antibody Name:** Anti GAPDH in vitro Antibody

**Description:** This monoclonal targets GAPDH

**Clone ID:** 6C5cc

**Antibody ID:** AB\_2858176

**Vendor:** Hytest

**Catalog Number:** 5G4cc-6C5cc

**Record Creation Time:** 20260218T234315+0000

**Record Last Update:** 20260226T000455+0000

---

### Ratings and Alerts

No rating or validation information has been found for Anti GAPDH in vitro Antibody.

No alerts have been found for Anti GAPDH in vitro Antibody.

## Data and Source Information

**Source:** [Antibody Registry](#)

## Usage and Citation Metrics

We found 6 mentions in open access literature.

**Listed below are recent publications.** The full list is available at [RRID](#) .

Ishino Y, et al. (2026) LPLAT11/MBOAT7-driven phosphatidylinositol remodeling ensures radial glial cell integrity in developing neocortex. *iScience*, 29(1), 114248.

Jäntti NZ, et al. (2025) EPLIN? controls integrin recycling from Rab21 endosomes to drive breast cancer cell migration. *Developmental cell*.

Kirkgöz K, et al. (2025) An intracellular recombinant single-chain variable antibody fragment as a new class of phosphodiesterase type 5 inhibitors. *British journal of pharmacology*.

Pieper NM, et al. (2024) Inhibition of bromodomain and extra-terminal proteins targets constitutively active NF?B and STAT signaling in lymphoma and influences the expression of the antiapoptotic proteins BCL2A1 and c-MYC. *Cell communication and signaling : CCS*, 22(1), 415.

De Jong KA, et al. (2023) Live-cell imaging identifies cAMP microdomains regulating ?-adrenoceptor-subtype-specific lipolytic responses in human white adipocytes. *Cell reports*, 42(5), 112433.

Peuhu E, et al. (2022) MYO10-filopodia support basement membranes at pre-invasive tumor boundaries. *Developmental cell*, 57(20), 2350.