

Resource Summary Report

Generated by [RRID](#) on Aug 6, 2026

Anti-dimethyl-Arginine, asymmetric (ASYM24)

RRID:AB_310596

Type: Antibody

Proper Citation

Millipore Cat# 07-414, RRID:AB_310596

Antibody Information

URL: http://antibodyregistry.org/AB_310596

Proper Citation: Millipore Cat# 07-414, RRID:AB_310596

Target Antigen: dimethyl-Arginine asymmetric (ASYM24)

Host Organism: rabbit

Clonality: polyclonal

Comments: seller recommendations: IgG; IgG WB; Western Blot

Antibody Name: Anti-dimethyl-Arginine, asymmetric (ASYM24)

Description: This polyclonal targets dimethyl-Arginine asymmetric (ASYM24)

Target Organism: h, m

Antibody ID: AB_310596

Vendor: Millipore

Catalog Number: 07-414

Record Creation Time: 20250820T064350+0000

Record Last Update: 20250821T004745+0000

Ratings and Alerts

No rating or validation information has been found for Anti-dimethyl-Arginine, asymmetric (ASYM24).

No alerts have been found for Anti-dimethyl-Arginine, asymmetric (ASYM24).

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We found 12 mentions in open access literature.

Listed below are recent publications. The full list is available at [RRID](#) .

Zhou L, et al. (2025) Mina53 catalyzes arginine demethylation of p53 to promote tumor growth. *Cell reports*, 44(2), 115242.

Gao J, et al. (2025) PIWI proteins tether the piRNA biogenesis machinery to mitochondria during mammalian spermatogenesis. *The EMBO journal*, 44(22), 6397.

Chang CH, et al. (2025) Nucleolar fibrillarin methyltransferase regulates systemic trafficking of a plant virus satellite RNA. *The Plant cell*, 37(10).

Gao J, et al. (2024) TDRD1 phase separation drives intermitochondrial cement assembly to promote piRNA biogenesis and fertility. *Developmental cell*, 59(20), 2704.

Wang L, et al. (2024) Arginine methylation-enabled FUS phase separation with SMN contributes to neuronal granule formation. *Cell reports*, 43(8), 114537.

Antal CE, et al. (2023) A super-enhancer-regulated RNA-binding protein cascade drives pancreatic cancer. *Nature communications*, 14(1), 5195.

Migazzi A, et al. (2021) Huntingtin-mediated axonal transport requires arginine methylation by PRMT6. *Cell reports*, 35(2), 108980.

Courchaine EM, et al. (2021) DMA-tudor interaction modules control the specificity of in vivo condensates. *Cell*, 184(14), 3612.

Huang T, et al. (2021) PRMT6 methylation of RCC1 regulates mitosis, tumorigenicity, and radiation response of glioblastoma stem cells. *Molecular cell*, 81(6), 1276.

Hu P, et al. (2019) Dual regulation of Arabidopsis AGO2 by arginine methylation. *Nature communications*, 10(1), 844.

Qamar S, et al. (2018) FUS Phase Separation Is Modulated by a Molecular Chaperone and Methylation of Arginine Cation-? Interactions. *Cell*, 173(3), 720.

Bao X, et al. (2017) CSNK1a1 Regulates PRMT1 to Maintain the Progenitor State in Self-Renewing Somatic Tissue. *Developmental cell*, 43(2), 227.