

Resource Summary Report

Generated by [RRID](#) on Sep 10, 2026

EGF antibody

RRID:AB_368486

Type: Antibody

Proper Citation

GeneTex Cat# GTX29695, RRID:AB_368486

Antibody Information

URL: http://antibodyregistry.org/AB_368486

Proper Citation: GeneTex Cat# GTX29695, RRID:AB_368486

Target Antigen: EGF antibody

Host Organism: rabbit

Clonality: polyclonal

Comments: Discontinued; manufacturer recommendations: IgG; IgG ELISA; Western Blot; Block/Neutralize/Inhibit; Neutralization: To yield one-half maximal inhibition [ND₅₀] of the biological activity of hEGF (2.0 ng/ml), a concentration of 0.25 - 0.35 µg/ml of this antibody is required. ELISA: To detect hEGF by direct ELISA (using 100µl/well antibody solution) a concentration of at least 0.5µg/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with compatible secondary reagents, allows the detection of 0.2 - 0.4 ng/well of recombinant hEGF. Western Blot: To detect hEGF by Western Blot analysis this antibody can be used at a concentration of 0.1 - 0.2 µg/ml. Used in conjunction with compatible secondary reagents the detection limit for recombinant hEGF is 1.5 - 3.0 ng/lane, under either

Antibody Name: EGF antibody

Description: This polyclonal targets EGF antibody

Target Organism: human

Antibody ID: AB_368486

Vendor: GeneTex

Catalog Number: GTX29695

Record Creation Time: 20231110T081302+0000

Ratings and Alerts

No rating or validation information has been found for EGF antibody.

Warning: Discontinued at GeneTex

Discontinued; manufacturer recommendations: IgG; IgG ELISA; Western Blot; Block/Neutralize/Inhibit; Neutralization: To yield one-half maximal inhibition [ND₅₀] of the biological activity of hEGF (2.0 ng/ml), a concentration of 0.25 - 0.35 µg/ml of this antibody is required. ELISA: To detect hEGF by direct ELISA (using 100µl/well antibody solution) a concentration of at least 0.5µg/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with compatible secondary reagents, allows the detection of 0.2 - 0.4 ng/well of recombinant hEGF. Western Blot: To detect hEGF by Western Blot analysis this antibody can be used at a concentration of 0.1 - 0.2 µg/ml. Used in conjunction with compatible secondary reagents the detection limit for recombinant hEGF is 1.5 - 3.0 ng/lane, under either

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.