

# Resource Summary Report

Generated by [RRID](#) on Sep 11, 2026

## E Cadherin antibody

RRID:AB\_368782

Type: Antibody

---

### Proper Citation

GeneTex Cat# GTX15148, RRID:AB\_368782

---

### Antibody Information

**URL:** [http://antibodyregistry.org/AB\\_368782](http://antibodyregistry.org/AB_368782)

**Proper Citation:** GeneTex Cat# GTX15148, RRID:AB\_368782

**Target Antigen:** E Cadherin antibody

**Host Organism:** rabbit

**Clonality:** polyclonal

**Comments:** Discontinued; manufacturer recommendations: IgG; IgG Immunohistochemistry - fixed; Immunohistochemistry; Immunohistochemistry, IHC (Formalin-fixed paraffin-embedded sections). IHC-P: Use at a dilution of 1/30 for 10 min at RT. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM citrate buffer, pH 6.0 for 10 min followed by cooling at RT for 20 min. Not tested in other applications. Optimal dilutions/concentrations should be determined by the end user.

**Antibody Name:** E Cadherin antibody

**Description:** This polyclonal targets E Cadherin antibody

**Antibody ID:** AB\_368782

**Vendor:** GeneTex

**Catalog Number:** GTX15148

**Record Creation Time:** 20231110T081236+0000

**Record Last Update:** 20260829T091412+0000

---

### Ratings and Alerts

No rating or validation information has been found for E Cadherin antibody.

**Warning:** Discontinued at GeneTex

Discontinued; manufacturer recommendations: IgG; IgG Immunohistochemistry - fixed; Immunohistochemistry; Immunohistochemistry, IHC (Formalin-fixed paraffin-embedded sections). IHC-P: Use at a dilution of 1/30 for 10 min at RT. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM citrate buffer, pH 6.0 for 10 min followed by cooling at RT for 20 min. Not tested in other applications. Optimal dilutions/concentrations should be determined by the end user.

---

## Data and Source Information

**Source:** [Antibody Registry](#)

---

## Usage and Citation Metrics

We have not found any literature mentions for this resource.