

Resource Summary Report

Generated by [RRID](#) on Aug 1, 2026

WDR5 antibody - ChIP Grade

RRID:AB_946146

Type: Antibody

Proper Citation

Abcam Cat# ab56919, RRID:AB_946146

Antibody Information

URL: http://antibodyregistry.org/AB_946146

Proper Citation: Abcam Cat# ab56919, RRID:AB_946146

Target Antigen: WDR5

Host Organism: mouse

Clonality: monoclonal

Comments: validation status unknown, seller recommendations provided in 2012: Western Blot; Western Blot

Antibody Name: WDR5 antibody - ChIP Grade

Description: This monoclonal targets WDR5

Target Organism: human

Antibody ID: AB_946146

Vendor: Abcam

Catalog Number: ab56919

Record Creation Time: 20250819T233335+0000

Record Last Update: 20250820T054243+0000

Ratings and Alerts

No rating or validation information has been found for WDR5 antibody - ChIP Grade.

No alerts have been found for WDR5 antibody - ChIP Grade.

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We found 6 mentions in open access literature.

Listed below are recent publications. The full list is available at [RRID](#) .

Yayli G, et al. (2023) ATAC and SAGA co-activator complexes utilize co-translational assembly, but their cellular localization properties and functions are distinct. Cell reports, 42(9), 113099.

Liao X, et al. (2023) Downregulation of BASP1 promotes temozolomide resistance in gliomas via epigenetic activation of the FBXO32/NF- κ B/MGMT axis. Molecular cancer research : MCR.

Fischer V, et al. (2021) The related coactivator complexes SAGA and ATAC control embryonic stem cell self-renewal through acetyltransferase-independent mechanisms. Cell reports, 36(8), 109598.

Luo H, et al. (2019) HOTTIP lncRNA Promotes Hematopoietic Stem Cell Self-Renewal Leading to AML-like Disease in Mice. Cancer cell, 36(6), 645.

Fang L, et al. (2018) SET1A-Mediated Mono-Methylation at K342 Regulates YAP Activation by Blocking Its Nuclear Export and Promotes Tumorigenesis. Cancer cell, 34(1), 103.

Hernandez C, et al. (2018) Dppa2/4 Facilitate Epigenetic Remodeling during Reprogramming to Pluripotency. Cell stem cell, 23(3), 396.