

Resource Summary Report

Generated by [RRID](#) on Jul 29, 2026

Anti-Rabbit IgG (H+L), made in goat

RRID:AB_2336197

Type: Antibody

Proper Citation

Vector Laboratories Cat# FI-1000, RRID:AB_2336197

Antibody Information

URL: http://antibodyregistry.org/AB_2336197

Proper Citation: Vector Laboratories Cat# FI-1000, RRID:AB_2336197

Target Antigen: IgG

Host Organism: goat

Clonality: polyclonal

Comments: Applications: IF, IHC, IHC-IF, ICC-IF, Flow cytometry/FACS, WB

Antibody Name: Anti-Rabbit IgG (H+L), made in goat

Description: This polyclonal targets IgG

Target Organism: rabbit

Antibody ID: AB_2336197

Vendor: Vector Laboratories

Catalog Number: FI-1000

Record Creation Time: 20250820T033949+0000

Record Last Update: 20250820T165104+0000

Ratings and Alerts

No rating or validation information has been found for Anti-Rabbit IgG (H+L), made in goat.

No alerts have been found for Anti-Rabbit IgG (H+L), made in goat.

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We found 27 mentions in open access literature.

Listed below are recent publications. The full list is available at [RRID](#) .

Shaimardanova AA, et al. (2024) Increasing α -hexosaminidase A activity using genetically modified mesenchymal stem cells. *Neural regeneration research*, 19(1), 212.

Mastroiacovo F, et al. (2022) Within the Ischemic Penumbra, Sub-Cellular Compartmentalization of Heat Shock Protein 70 Overlaps with Autophagy Proteins and Fails to Merge with Lysosomes. *Molecules (Basel, Switzerland)*, 27(10).

Wu S, et al. (2021) Endoplasmic reticulum associated degradation is required for maintaining endoplasmic reticulum homeostasis and viability of mature Schwann cells in adults. *Glia*, 69(2), 489.

Flowers BM, et al. (2021) Cell of Origin Influences Pancreatic Cancer Subtype. *Cancer discovery*, 11(3), 660.

Ishikita A, et al. (2021) GFAT2 mediates cardiac hypertrophy through HBP-O-GlcNAcylation-Akt pathway. *iScience*, 24(12), 103517.

Nagayach A, et al. (2021) Connected neurons in multiple neocortical areas, comprising parallel circuits, encode essential information for visual shape learning. *Journal of chemical neuroanatomy*, 118, 102024.

Katarzyna Greda A, et al. (2021) Hyaluronidase inhibition accelerates functional recovery from stroke in the mouse brain. *Journal of neurochemistry*, 157(3), 781.

Wang TE, et al. (2021) Secretory mouse quiescin sulfhydryl oxidase 1 aggregates defected human and mouse spermatozoa in vitro and in vivo. *iScience*, 24(10), 103167.

Stone S, et al. (2020) The UPR preserves mature oligodendrocyte viability and function in adults by regulating autophagy of PLP. *JCI insight*, 5(5).

Cheng Y, et al. (2020) Neuroprotective actions of leptin facilitated through balancing mitochondrial morphology and improving mitochondrial function. *Journal of neurochemistry*, 155(2), 191.

Tang X, et al. (2020) Using inducible lentiviral vectors to generate induced pluripotent stem cell line ZOCi001-A from peripheral blood cells of a patient with CRB1^{-/-} retinitis pigmentosa. *Stem cell research*, 45, 101817.

Tang X, et al. (2020) Generation of a homozygous CRISPR/Cas9-mediated knockout H9 hESC subline for the CRB1 locus. *Stem cell research*, 49, 102057.

Park H, et al. (2020) Generation of a human induced pluripotent stem cell line from a patient with Leber congenital amaurosis. Stem cell research, 43, 101725.

Bieging-Rolett KT, et al. (2020) Zmat3 Is a Key Splicing Regulator in the p53 Tumor Suppression Program. Molecular cell, 80(3), 452.

Lei Z, et al. (2020) NF- κ B Activation Accounts for the Cytoprotective Effects of PERK Activation on Oligodendrocytes during EAE. The Journal of neuroscience : the official journal of the Society for Neuroscience, 40(33), 6444.

Wu S, et al. (2020) The Integrated UPR and ERAD in Oligodendrocytes Maintain Myelin Thickness in Adults by Regulating Myelin Protein Translation. The Journal of neuroscience : the official journal of the Society for Neuroscience, 40(43), 8214.

Bowen ME, et al. (2019) The Spatiotemporal Pattern and Intensity of p53 Activation Dictates Phenotypic Diversity in p53-Driven Developmental Syndromes. Developmental cell, 50(2), 212.

Kim BY, et al. (2019) Generation of a patient-specific induced pluripotent stem cell line, KSCBi006-A, for osteogenesis imperfecta type I with the COL1A1, c.3162delT mutation. Stem cell research, 41, 101622.

Park H, et al. (2019) Generation of human induced pluripotent stem cells from peripheral blood mononuclear cells of a Senior-Loken syndrome patient. Stem cell research, 41, 101648.

Nagayach A, et al. (2019) Separate Gene Transfers into Pre- and Postsynaptic Neocortical Neurons Connected by mGluR5-Containing Synapses. Journal of molecular neuroscience : MN, 68(4), 549.