

Resource Summary Report

Generated by [RRID](#) on Aug 10, 2026

CD133/1 (AC133)-APC, human

RRID:AB_2660883

Type: Antibody

Proper Citation

Miltenyi Biotec Cat# 130-098-829, RRID:AB_2660883

Antibody Information

URL: http://antibodyregistry.org/AB_2660883

Proper Citation: Miltenyi Biotec Cat# 130-098-829, RRID:AB_2660883

Target Antigen: CD133/1 (AC133)

Host Organism: mouse

Clonality: monoclonal

Comments: Discontinued: 6-2018; Target Distribution brain, endothelial cells, epithelial cells, heart, hematopoietic stem cells, kidney, liver, lung, pancreas, placenta, ES and iPS cells, red blood cells; target type CD markers; tested applications MACS Flow Cytometry; quantity:

Info: This product is discontinued and reformatted to a higher concentration for optimized use in multicolor flow cytometry panels. The replacement product cat # is 130-113-668. (RRID:AB_2726210).

Antibody Name: CD133/1 (AC133)-APC, human

Description: This monoclonal targets CD133/1 (AC133)

Target Organism: human

Clone ID: AC133

Antibody ID: AB_2660883

Vendor: Miltenyi Biotec

Catalog Number: 130-098-829

Record Creation Time: 20250820T053457+0000

Record Last Update: 20250820T215551+0000

Ratings and Alerts

No rating or validation information has been found for CD133/1 (AC133)-APC, human.

Warning: Discontinued: 2021

Discontinued: 6-2018; Target Distribution brain, endothelial cells, epithelial cells, heart, hematopoietic stem cells, kidney, liver, lung, pancreas, placenta, ES and iPS cells, red blood cells; target type CD markers; tested applications MACS Flow Cytometry; quantity: Info: This product is discontinued and reformatted to a higher concentration for optimized use in multicolor flow cytometry panels. The replacement product cat # is 130-113-668. (RRID:AB_2726210).

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We found 3 mentions in open access literature.

Listed below are recent publications. The full list is available at [RRID](#) .

Tartaglia JT, et al. (2024) Mobilization of Endogenous CD34+/CD133+ Endothelial Progenitor Cells by Enhanced External Counter Pulsation for Treatment of Refractory Angina. International journal of molecular sciences, 25(18).

Hallett JM, et al. (2022) Human biliary epithelial cells from discarded donor livers rescue bile duct structure and function in a mouse model of biliary disease. Cell stem cell, 29(3), 355.

Liau BB, et al. (2017) Adaptive Chromatin Remodeling Drives Glioblastoma Stem Cell Plasticity and Drug Tolerance. Cell stem cell, 20(2), 233.