

Resource Summary Report

Generated by [RRID](#) on Sep 10, 2026

p27 KIP1 [DCS-72.F6] antibody

RRID:AB_369307

Type: Antibody

Proper Citation

GeneTex Cat# GTX23928, RRID:AB_369307

Antibody Information

URL: http://antibodyregistry.org/AB_369307

Proper Citation: GeneTex Cat# GTX23928, RRID:AB_369307

Target Antigen: p27 KIP1 [DCS-72.F6] antibody

Host Organism: mouse

Clonality: monoclonal

Comments: Discontinued; manufacturer recommendations: Immunohistochemistry - frozen; Immunoprecipitation; Immunohistochemistry; Western Blot; Immunohistochemistry - fixed; Immunoprecipitation: use at 1/50. Western Blotting: use at 1/50 to 1/100. This antibody may be diluted to a titer of 1/25-1/50 in an ABC method. However, depending on the staining conditions employed, we suggest that final dilution should be determined by the user. We suggest an incubation period of 30 minutes at room temperature. Formalin fixed paraffin embedded tissue sections require high temperature antigen unmasking with 10 mM citrate buffer, pH 6.0 prior to immunostaining. employed, we suggest that final dilution should be determined by the user. We suggest an incubation period of 30 minutes at room temperature. Formalin fixed paraffin embedded tissue sections require h

Antibody Name: p27 KIP1 [DCS-72.F6] antibody

Description: This monoclonal targets p27 KIP1 [DCS-72.F6] antibody

Target Organism: mouse, human

Antibody ID: AB_369307

Vendor: GeneTex

Catalog Number: GTX23928

Record Creation Time: 20231110T081302+0000

Ratings and Alerts

No rating or validation information has been found for p27 KIP1 [DCS-72.F6] antibody.

Warning: Discontinued at GeneTex

Discontinued; manufacturer recommendations: Immunohistochemistry - frozen; Immunoprecipitation; Immunohistochemistry; Western Blot; Immunohistochemistry - fixed; Immunoprecipitation: use at 1/50. Western Blotting: use at 1/50 to 1/100. This antibody may be diluted to a titer of 1/25-1/50 in an ABC method. However, depending on the staining conditions employed, we suggest that final dilution should be determined by the user. We suggest an incubation period of 30 minutes at room temperature. Formalin fixed paraffin embedded tissue sections require high temperature antigen unmasking with 10 mM citrate buffer, pH 6.0 prior to immunostaining. employed, we suggest that final dilution should be determined by the user. We suggest an incubation period of 30 minutes at room temperature. Formalin fixed paraffin embedded tissue sections require h

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.