

# Resource Summary Report

Generated by [RRID](#) on Aug 5, 2026

## Cytokeratin 18 [CY90] antibody (ready-to-use)

RRID:AB\_376824

Type: Antibody

### Proper Citation

GeneTex Cat# GTX75467, RRID:AB\_376824

### Antibody Information

**URL:** [http://antibodyregistry.org/AB\\_376824](http://antibodyregistry.org/AB_376824)

**Proper Citation:** GeneTex Cat# GTX75467, RRID:AB\_376824

**Target Antigen:** Cytokeratin 18 [CY90] antibody (ready-to-use)

**Host Organism:** mouse

**Clonality:** monoclonal

**Comments:** Discontinued; manufacturer recommendations: Immunohistochemistry - fixed; Immunohistochemistry; Immunohistochemistry - frozen; IHC (Frozen sections), IHC (Formalin-fixed paraffin-embedded sections). The usefulness of this product in other applications has not been determined., IHC-Fr, IHC-P, For IHC-Fr and IHC-P: Use at a dilution of neat - 1:10. Not tested in other applications. Optimal concentrations/dilutions should be determined by the researcher. \*This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Serotec recommend STUF (HIS003) for this procedure.

**Antibody Name:** Cytokeratin 18 [CY90] antibody (ready-to-use)

**Description:** This monoclonal targets Cytokeratin 18 [CY90] antibody (ready-to-use)

**Target Organism:** human

**Antibody ID:** AB\_376824

**Vendor:** GeneTex

**Catalog Number:** GTX75467

**Record Creation Time:** 20250820T011354+0000

**Record Last Update:** 20250820T102240+0000

---

## Ratings and Alerts

No rating or validation information has been found for Cytokeratin 18 [CY90] antibody (ready-to-use).

**Warning:** Discontinued at GeneTex

Discontinued; manufacturer recommendations: Immunohistochemistry - fixed; Immunohistochemistry; Immunohistochemistry - frozen; IHC (Frozen sections), IHC (Formalin-fixed paraffin-embedded sections). The usefulness of this product in other applications has not been determined., IHC-Fr, IHC-P, For IHC-Fr and IHC-P: Use at a dilution of neat - 1:10. Not tested in other applications. Optimal concentrations/dilutions should be determined by the researcher. \*This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Serotec recommend STUF (HIS003) for this procedure.

---

## Data and Source Information

**Source:** [Antibody Registry](#)

---

## Usage and Citation Metrics

We have not found any literature mentions for this resource.