

Resource Summary Report

Generated by [RRID](#) on Aug 30, 2026

GFP-C1-PLCdelta-PH

RRID:Addgene_21179

Type: Plasmid

Proper Citation

RRID:Addgene_21179

Plasmid Information

URL: <http://www.addgene.org/21179>

Proper Citation: RRID:Addgene_21179

Insert Name: PLCdelta

Organism: Homo sapiens

Bacterial Resistance: Kanamycin

Defining Citation: [PMID:9512420](#)

Vector Backbone Description: Backbone Marker:Clontech; Backbone Size:4700; Vector Backbone:C1; Vector Types:Mammalian Expression; Bacterial Resistance:Kanamycin

Comments: This construct encodes a PI(4,5)P2 lipid selective PH domain that can be used as a fluorescent translocation biosensor to monitor changes or local differences in the concentration of plasma membrane PI(4,5)P2 lipids (Stauffer TP, Ahn S, Meyer T; Curr Biol. 1998 8:343-6). It is best used in combination with a second fluorescent protein that uniformly marks the plasma membrane (the last 20 amino acids "CAAX" region of KRas is often used as a PM marker). Relative changes in cell-wide or local plasma membrane PI(4,5)P2 can be determined by using a mask for the PM region and by dividing the two background subtracted PH domain and CAAX images. At high expression levels, this construct interferes with PI(4,5)P2 signaling and alters cortical actin (Raucher D et al Cell. 2000 100:221-8).

Plasmid Name: GFP-C1-PLCdelta-PH

Relevant Mutation: PH domain of human PLC delta, amino acids Met1–Ile175. (Lab plasmid WC 0003)

Record Creation Time: 20220422T222054+0000

Record Last Update: 20260829T080202+0000

Ratings and Alerts

No rating or validation information has been found for GFP-C1-PLCdelta-PH.

No alerts have been found for GFP-C1-PLCdelta-PH.

Data and Source Information

Source: [Addgene](#)

Usage and Citation Metrics

We found 96 mentions in open access literature.

Listed below are recent publications. The full list is available at [RRID](#) .

Bott CJ, et al. (2026) Lipid transfer protein ORP3 mediates lysosomal repair via LC3B and ubiquitin-TAK1-p38 signaling. bioRxiv : the preprint server for biology.

Ferreira JV, et al. (2026) LAMP2A regulates endosomal protein composition and membrane identity in exosome biogenesis. iScience, 29(1), 114305.

Sejd JR, et al. (2026) Chemical Genetic Screen Identifies PSD3 as a Direct Substrate of NUA1 that Regulates Dendritic Spine Maturation. bioRxiv : the preprint server for biology.

Abdelnaser RA, et al. (2025) Identification of TNFAIP2 as a unique cellular regulator of CSF-1 receptor activation. Life science alliance, 8(5).

Verpoort B, et al. (2025) A postsynaptic GPR158-PLCXD2 complex controls spine apparatus abundance and dendritic spine maturation. Developmental cell.

Taciak B, et al. (2025) Harnessing macrophage-drug conjugates for allogeneic cell-based therapy of solid tumors via the TRAIN mechanism. Nature communications, 16(1), 1327.

Jha A, et al. (2025) CD44 and Ezrin restrict EGF receptor mobility to generate a novel spatial arrangement of cytoskeletal signaling modules driving bleb-based migration. bioRxiv : the preprint server for biology.

Ozaki-Noma R, et al. (2025) Positive-Type Reversibly Photoswitching Red Fluorescent Protein for Dual-Color Superresolution Imaging with Single Light Exposure for Off-Switching. ACS nano, 19(7), 7188.

Reinders M, et al. (2025) ATXN3 regulates lysosome regeneration after damage by targeting K48-K63-branched ubiquitin chains. The EMBO journal, 44(18), 5086.

Christodoulou N, et al. (2025) Calcium transients regulate the apical emergence of basally located progenitors during Xenopus skin development. Nature communications, 16(1), 6650.

Hiyoshi M, et al. (2025) M-Sec promotes the accumulation of intracellular HTLV-1 Gag puncta and the incorporation of Env into viral particles. PLoS pathogens, 21(1), e1012919.

Elsayyid MW, et al. (2025) Phosphatidylinositol 4,5-bisphosphate impacts extracellular vesicle shedding from C. elegans ciliated sensory neurons. Journal of cell science, 138(20).

Elsayyid M, et al. (2025) Phosphatidylinositol (4,5)-bisphosphate Impacts Ectosome Shedding from C. elegans Ciliated Sensory Neurons. bioRxiv : the preprint server for biology.

Dai A, et al. (2025) Multiple interactions mediate the localization of BLTP2 at ER-PM contacts to control plasma membrane dynamics. bioRxiv : the preprint server for biology.

He F, et al. (2025) Roles of class III phosphatidylinositol 3-kinase, Vps34, in phagocytosis, autophagy, and endocytosis in retinal pigmented epithelium. iScience, 28(5), 112371.

Lee C, et al. (2025) Optogenetic storage and release of protein and mRNA in live cells and animals. Nature communications, 16(1), 6230.

Ajiki M, et al. (2024) ORP9-PH domain-based fluorescent reporters for visualizing phosphatidylinositol 4-phosphate dynamics in living cells. RSC chemical biology, 5(6), 544.

George K, et al. (2024) Robust GRK2/3/6-dependent desensitization of oxytocin receptor in neurons. iScience, 27(6), 110047.

Di Meo D, et al. (2024) Pip5k1 γ regulates axon formation by limiting Rap1 activity. Life science alliance, 7(5).

Brands J, et al. (2024) A molecular mechanism to diversify Ca²⁺ signaling downstream of Gs protein-coupled receptors. Nature communications, 15(1), 7684.