

# Resource Summary Report

Generated by [RRID](#) on Aug 11, 2026

## Thermo Fisher: TSQ 9000 Triple Quadrupole GC-MS/MS system

RRID:SCR\_026311

Type: Tool

### Proper Citation

Thermo Fisher: TSQ 9000 Triple Quadrupole GC-MS/MS system (RRID:SCR\_026311)

### Resource Information

**URL:** <https://www.thermofisher.com/order/catalog/product/TTS90BOP>

**Proper Citation:** Thermo Fisher: TSQ 9000 Triple Quadrupole GC-MS/MS system (RRID:SCR\_026311)

**Description:** Triple quadrupole mass spectrometry system to study liquid, gaseous or solid samples. Mass analyzer enables absolute quantification of low-concentration metabolites. Designed for complex sample matrices in high-throughput applications. ExtractaBrite or AEI configuration are used to maximize uptime between routine maintenance interventions. The wireless ExtractaBrite ion source is heated throughout and includes the first RF region, commonly the first place susceptible to ion burn. This design ensures high matrix tolerance for high throughput applications. The AEI source is also heated throughout, with RF lenses and benefits from an innate ultra-robustness due to a more tightly focused ion beam, limiting ion burn on the source lenses.

**Resource Type:** instrument resource

**Keywords:** Triple quadrupole mass spectrometry, liquid, gaseous, solid, sample, mass analyzer, quantification of low-concentration metabolites,

**Specification URL:** [https://raw.githubusercontent.com/SciCrunch/RRID-Instruments/main/PDF/SCR\\_026311.pdf](https://raw.githubusercontent.com/SciCrunch/RRID-Instruments/main/PDF/SCR_026311.pdf)

**Resource Name:** Thermo Fisher: TSQ 9000 Triple Quadrupole GC-MS/MS system

**Resource ID:** SCR\_026311

**Alternate IDs:** Model\_Number\_TSQ\_9000

**Alternate URLs:**

<https://www.thermofisher.cn/order/catalog/product/TTS90MTS?SID=srch-srp-TTS90MTS>

**Record Creation Time:** 20250118T053338+0000

**Record Last Update:** 20260808T190738+0000

---

## Ratings and Alerts

No rating or validation information has been found for Thermo Fisher: TSQ 9000 Triple Quadrupole GC-MS/MS system.

No alerts have been found for Thermo Fisher: TSQ 9000 Triple Quadrupole GC-MS/MS system.

---

## Data and Source Information

**Source:** [SciCrunch Registry](#)

---

## Usage and Citation Metrics

We have not found any literature mentions for this resource.